




ENMED® ENM - 410 Disposable Nonwoven Bouffant Cap

The ENMED ENM - 410 Disposable Nonwoven Bouffant Cap offers Adjustable soft thread elastic for comfort fittings, sonic sealing, attractive look, avoids escape of hair, dandruff and micro organisms during process ENM-410 is made from spun bond polymer non-woven fabric, spun lace fabrics as the main raw material. This material is preferred for its light weight. Its gravity is lighter than Hydrophobic PP fibers so it has excellent water repellent and air permeability resulting in providing comfort & dryness resulting in anti bacterial and anti fungal resistance.

INTENDED USE

Non-sterile bouffant caps intended to be worn by health care personnel during various procedures to protect patient/materials and the operating personnel from transfer of microorganisms, particulates, etc.



Color and Size Table

Color	Size
	46 / 52 cm (19 / 21 Inch)
Blue	A8008-PLAS










Dimension Table

Description		Dimension (cm)
		One size
B	Shoe Cover Height	15
A	Shoe Cover Natural Length	38
C	Opening Stretched Length	37
C	Opening Natural Length	9

Storage and Disposal

- Store in dry, clean conditions in original packaging, away from direct sunlight, sources of high temperature, and solvent vapours
- Store within the temperature range -5°C to +25°C and with relative humidity below 80%
- Expected shelf life is three years from date of manufacture when stored as stated
- Replace garments if damaged, heavily contaminated or in accordance with local work practice or regulations
- Handle and dispose of contaminated garments with care and in accordance with applicable regulations

Limited Use

	DO NOT RE USE		DO NOT BLEACH		DO NOT IRON
	DO NOT DRY CLEAN		FLAMMABLE KEEP AWAY FROM SPARKS OF FLAMES		CONSULT INSTRUCTIONS FOR USE
	DO NOT WASH LAUNDERING IMPACTS UPON PROTECTIVE PERFORMANCE		DO NOT MACHINE DRY		CAUTION